# GEMINUS<sup>®</sup> Volar Distal Radius Plating System

## INSTRUCTIONS FOR USE

**R:** For use by physicians only. Federal Law restricts this device to sale by or on the order of a physician.

### Failure to follow instructions may lead to patient injury.

This package insert is designed to provide Instructions for Use of the GEMINUS Volar Plating System; it is not a reference to surgical techniques.

### SYMBOLS

|   |                                 |   |   |
|---|---------------------------------|---|---|
| <b>MATL:</b>  | MATERIAL                        | <b>CoCr:</b>  | COBALT CHROMIUM ALLOY                               |
| <b>MADE IN:</b>   | MADE IN <<COUNTRY>>             | <b>Ti:</b>  | TITANIUM ALLOY                                      |
| <b>QTY:</b>   | QUANTITY                        | <b>SS, SST:</b>   | STAINLESS STEEL                                     |
|    | DO NOT REUSE (SINGLE USE)       |    | CAUTION or ATTENTION, SEE INSTRUCTIONS FOR USE      |
|    | USE BY (EXPIRATION DATE)        |    | CONSULT INSTRUCTIONS FOR USE                        |
|    | BATCH CODE                      |    | MANUFACTURER  |
|    | STERILIZED USING ETHYLENE OXIDE |    | TEMPERATURE LIMITATION                              |
|    | STERILIZED USING IRRADIATION    |    | AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY |
|   | NON STERILE PRODUCT             |  | DO NOT USE IF PACKAGE IS DAMAGED                    |
|  | CATALOG NUMBER                  |   |   |

### Description:

The Skeletal Dynamics GEMINUS Volar Plating System contains bone plates for the repair of distal volar radial fractures. Included in the set are titanium bone screws, fixation pegs, fragment plates, and specialized instrumentation. Also included are a Hook Plate Extension to buttress a volar marginal fragment, and cannulated cobalt chrome polyaxial locking screws for trajectories different than those of the fixed angled bone plates.

The GEMINUS Volar Plates are available in various sizes and are made of medical grade titanium alloy. Cortical screws affix the plate to the diaphysis and fixed angle pegs are used for distal bone fragments. The system is provided non-sterile and is sterilized in the user facility.

The GEMINUS Volar Plating System is comprised of:

- Titanium alloy plates, washers and screws
- CoCr Cannulated Polyaxial Locking Screw (PLS)
- Stainless steel K-wires (for provisional fixation; not for implantation)
- System specific instrumentation

### Indications:

The GEMINUS Volar Plating System is intended for the fixation of fractures and osteotomies involving the volar distal radius.

### Contraindications:

Prior to using the GEMINUS Volar Plating System, ensure that none of the following patient conditions are present: active or latent infection, sepsis, insufficient quantity or quality of bone and/or soft tissue, material sensitivity, or patients who are unwilling or incapable of following post operative care instructions.

### **⚠ Warnings:**

- All screws must be implanted and fully tightened into the plate to maintain the integrity and strength of the finished construct. If the screws are not attached and/or fully tightened, a non-union, delayed union or construct failure may occur.
- The information in this document should be shared with the patient.
- The patient should be informed about the importance of following the post operative rehabilitation prescribed in order to fully understand the possible limitations in activities of daily living. The patient must be warned that failure to follow postoperative care instructions may cause the implant or treatment to fail.
- Potential GEMINUS<sup>®</sup> Volar Plating System construct failures such as stress fractures of the bones, loosening of the construct and/or fixation, delayed fusion, non-fusion, or incomplete healing may occur as a result of non compliance to post operative rehabilitation, excessive wrist activities or construct overloading.
- DO NOT reuse any of the GEMINUS Volar Plating System implantable components. Reuse may compromise the structural integrity of the construct and/or lead to failure or infection, which may result in patient injury.
- DO NOT open the volar capsule as it may devascularize fracture fragments and destabilize the volar wrist ligaments.
- Use only one 2.7mm Peg (High Compression or Fully Threaded, Non Locking) in each head of the GEMINUS Volar Plate.
- Use only one 2.5mm PLS in each head of the GEMINUS Volar Plate.
- DO NOT use the PLS in the most distal hole(s) on the lunate head of the GEMINUS Volar Plate.
- GEMINUS Drill Blocks are only compatible with GEMINUS Volar Plates containing a Gold PDG in the shaft.

### **⚠ Precautions:**

- Protect the GEMINUS Volar Plating System's implantable components against scratching or nicking. Such stress concentration can lead to implant failure.
- Before using the GEMINUS Volar Plating System, inspect all implants and instruments for wear, disfiguration and physical damage. If evidence of wear, disfiguration or physical damage is found, DO NOT use and contact your local Skeletal Dynamics representative or the Skeletal Dynamics Customer Care Department.
- DO NOT permanently implant the Skeletal Dynamics K-Wires; they are intended to be used during provisional fixation of the GEMINUS Volar Plate.
- DO NOT permanently implant the pre-loaded Drill Guides, Drill Blocks, or A.I.M.ing Guides; they are intended to be removed prior to peg insertion.
- The GEMINUS Volar Plating System has not been evaluated for safety and compatibility in the MR environment; nor has it been tested for heating or migration in the MR environment.
- DO NOT use peg/screw lengths that will excessively protrude through the far cortex as it may result in soft tissue irritation.
- The maximum angulation of the PLS should not exceed 10° from the trajectory of the respective hole.
- The Non-locking Threaded Pegs are NOT intended to provide subchondral support. Their use should be limited to capture remote bone fragments where partially or fully threaded pegs cannot be used.
- The Skeletal Dynamics GEMINUS Volar Plating System is to be used only with Skeletal Dynamics instruments, implants and accessories.
- Dispose of contaminated implants and instruments per established facility guidelines and protocols.
- Accuracy of Depth, Gap and Screw Gauges are within  $\pm 0.25$ mm.
- Caution should be taken for interference to pacemakers during electrocautery or by uncertified drills.
- Seek medical help immediately if implant malfunctions.
- To maintain traceability of the GEMINUS Volar Plating System implantable components, you must record each of the respective components LOT numbers into the patient records post implantation.

### **Potential Adverse Events:**

The following are potential risks that have been associated with wrist surgery: infection, nonunion, persistent pain, stiffness of the fingers, loosening or migration of the implants resulting in mal-alignment.

### **Directions for Use:**

The GEMINUS Volar Plating System should only be used by surgeons who have experience with this system. Each surgeon must evaluate the appropriateness for the use of the GEMINUS Volar Plating System based on their clinical experiences.

Please refer to the GEMINUS Volar Plating System's Surgical Technique Guide to review the surgical approach as described by Jorge L. Orbay, M.D. of the *Miami Hand and Upper Extremity Institute* located in Miami, Florida.

## Cleaning:

The GEMINUS® Volar Plating System instrumentation must be cleaned to achieve sterilization. The recommended manual cleaning instructions are set forth below. The user must validate other cleaning methods.

1. Disassemble instrumentation, if applicable.
2. Rinse components thoroughly under running cool tap water. While rinsing, use a soft bristle brush to loosen and remove as much visible soil as possible from components.
3. Soak components in a neutral enzymatic cleaner for a minimum of ten (10) minutes. Components must be fully immersed in the cleaner. Follow the cleaner manufacturer's instructions for cleaner preparation and exposure time.
4. Thoroughly rinse the components with cool water. While rinsing, use soft bristle brushes, pipettes or a water jet to clean out lumens, holes, and other challenging features.
5. Manually scrub the components thoroughly in newly made, clean, neutral pH enzymatic cleaner using soft bristle brushes or pipettes. All lumens, holes, hinged components, mating surfaces, and crevices, and challenging components should be thoroughly scrubbed. Actuate all moveable features and expose all areas to cleaner and to the brush or pipette.
6. Rinse components thoroughly with deionized or purified water; using pipettes or a water jet to clean out lumens, holes, and other hard to reach or challenging features. Actuate all movable features to fully irrigate all areas.
7. Visually inspect components for soil. Repeat the cleaning procedure until no visible soil remains on the components.
8. Perform a final rinse on the components using deionized water or purified water.
9. Dry the clean components using compressed air or a soft, lint free, clean cloth.

## Functional Checks should be performed where possible:

1. Mating devices should be checked for proper assembly.
2. Reusable devices with moving parts should be operated to check correct operation (medical grade lubricant suitable for steam sterilization can be applied as required).
3. Rotating instruments (e.g. drill bits, reamers) should be checked for straightness. This can be achieved by rolling the instrument on a flat surface.

**Note:** The useful life of these devices is dependent on many factors including, but not limited to the method and duration of each use and the handling of the devices between uses. Routine and careful inspection and functional testing of the device is the best method of determining the serviceable life span for the medical device.

## Sterilization:

The Skeletal Dynamics GEMINUS Volar Plating System is provided non sterile. This system is intended for steam sterilization at the healthcare facility.

1. Place all components and accessories into the designated areas of the sterilization tray
2. Steam sterilization may be accomplished using one of the cycles shown below:

| Cycle Type           | Temperature   | Duration             | Drying Time |
|----------------------|---------------|----------------------|-------------|
| Pre-Vacuum Autoclave | 270°F (132°C) | 4 minutes (wrapped)  | 20 minutes  |
| Gravity Autoclave    | 270°F (132°C) | 15 minutes (wrapped) | 20 minutes  |

- Follow ANSI/AAMI ST79:2006 - Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to the requirements of ANSI/AAMI ST79:2006 - Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Usage of an FDA approved wrap or sterilization container is required.
- Subsequent instrument sterilization needs to be performed in the tray system provided. For reuse and sterilization, instruments should be arranged within the tray system in the manner supplied by the company.

## Handling and Storage:

When not in use, store the clean and disinfected GEMINUS Volar Plating System within the Sterilization Tray. Prior to use, inspect the instrumentation for serviceability.

## Disclaimer of Warranty and Limited Remedies:

Skeletal Dynamics, LLC makes no express or implied warranty, including any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described in this publication. Skeletal Dynamics, LLC shall not be liable under any circumstances for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has authority to bind Skeletal Dynamics, LLC to any representation or warranty except as specifically set forth in this publication. Descriptions or specifications provided by Skeletal Dynamics, LLC in any publication are only included to generally describe the product when manufactured and do not constitute any express warranties.

**GEMINUS® Volar Plating System Ordering Information: GMN-FSP-SYS****Catalog #****Nomenclature**

| <b>Volar Distal Radius Plates (Ti)</b> |  |
|--|--|
| GMN-RTN-3HL                            | GEMINUS Volar Distal Radius Plate, Narrow, 3 Hole, Right   |
| GMN-LTN-3HL                            | GEMINUS Volar Distal Radius Plate, Narrow, 3 Hole, Left    |
| GMN-RTN-4HL                            | GEMINUS Volar Distal Radius Plate, Narrow, 4 Hole, Right   |
| GMN-LTN-4HL                            | GEMINUS Volar Distal Radius Plate, Narrow, 4 Hole, Left    |
| GMN-RTS-3HL                            | GEMINUS Volar Distal Radius Plate, Standard, 3 Hole, Right |
| GMN-LTS-3HL                            | GEMINUS Volar Distal Radius Plate, Standard, 3 Hole, Left  |
| GMN-RTS-4HL                            | GEMINUS Volar Distal Radius Plate, Standard, 4 Hole, Right |
| GMN-LTS-4HL                            | GEMINUS Volar Distal Radius Plate, Standard, 4 Hole, Left  |
| GMN-RTS-7HL                            | GEMINUS Volar Distal Radius Plate, Standard, 7 Hole, Right |
| GMN-LTS-7HL                            | GEMINUS Volar Distal Radius Plate, Standard, 7 Hole, Left  |
| GMN-RTW-4HL                            | GEMINUS Volar Distal Radius Plate, Wide, 4 Hole, Right     |
| GMN-LTW-4HL                            | GEMINUS Volar Distal Radius Plate, Wide, 4 Hole, Left      |
| GMN-HP                                 | GEMINUS Hook Plate   |
| GMN-HP-SCRW                            | GEMINUS Hook Plate, Screw                                  |
| WBTN-2750-T                            | Washer, Button, Inside 2.7mm x Outside 5.0mm, Ti           |
| <b>Smooth Pegs, Locking (Ti)</b>       |  |
| SPLS-20100-TS                          | Smooth Peg, Locking, 2.0mm x 10mm                          |
| SPLS-20120-TS                          | Smooth Peg, Locking, 2.0mm x 12mm                          |
| SPLS-20140-TS                          | Smooth Peg, Locking, 2.0mm x 14mm                          |
| SPLS-20160-TS                          | Smooth Peg, Locking, 2.0mm x 16mm                          |
| SPLS-20170-TS                          | Smooth Peg, Locking, 2.0mm x 17mm                          |
| SPLS-20180-TS                          | Smooth Peg, Locking, 2.0mm x 18mm                          |
| SPLS-20190-TS                          | Smooth Peg, Locking, 2.0mm x 19mm                          |
| SPLS-20200-TS                          | Smooth Peg, Locking, 2.0mm x 20mm                          |
| SPLS-20210-TS                          | Smooth Peg, Locking, 2.0mm x 21mm                          |
| SPLS-20220-TS                          | Smooth Peg, Locking, 2.0mm x 22mm                          |
| SPLS-20230-TS                          | Smooth Peg, Locking, 2.0mm x 23mm                          |
| SPLS-20240-TS                          | Smooth Peg, Locking, 2.0mm x 24mm                          |
| SPLS-20260-TS                          | Smooth Peg, Locking, 2.0mm x 26mm                          |
| SPLS-20280-TS                          | Smooth Peg, Locking, 2.0mm x 28mm                          |
| SPLS-20300-TS                          | Smooth Peg, Locking, 2.0mm x 30mm                          |
| SPLS-20320-TS                          | Smooth Peg, Locking, 2.0mm x 32mm                          |
| <b>Threaded Peg, Locking (Ti)</b>      |  |
| TPLS-23100-TS                          | Threaded Peg, Locking, 2.3mm x 10mm                        |
| TPLS-23120-TS                          | Threaded Peg, Locking, 2.3mm x 12mm                        |
| TPLS-23140-TS                          | Threaded Peg, Locking, 2.3mm x 14mm                        |
| TPLS-23160-TS                          | Threaded Peg, Locking, 2.3mm x 16mm                        |
| TPLS-23170-TS                          | Threaded Peg, Locking, 2.3mm x 17mm                        |
| TPLS-23180-TS                          | Threaded Peg, Locking, 2.3mm x 18mm                        |
| TPLS-23190-TS                          | Threaded Peg, Locking, 2.3mm x 19mm                        |
| TPLS-23200-TS                          | Threaded Peg, Locking, 2.3mm x 20mm                        |
| TPLS-23210-TS                          | Threaded Peg, Locking, 2.3mm x 21mm                        |
| TPLS-23220-TS                          | Threaded Peg, Locking, 2.3mm x 22mm                        |
| TPLS-23230-TS                          | Threaded Peg, Locking, 2.3mm x 23mm                        |
| TPLS-23240-TS                          | Threaded Peg, Locking, 2.3mm x 24mm                        |
| TPLS-23260-TS                          | Threaded Peg, Locking, 2.3mm x 26mm                        |
| TPLS-23280-TS                          | Threaded Peg, Locking, 2.3mm x 28mm                        |
| TPLS-23300-TS                          | Threaded Peg, Locking, 2.3mm x 30mm                        |
| TPLS-23320-TS                          | Threaded Peg, Locking, 2.3mm x 32mm                        |

| <b>High Compression Pegs, Locking (Ti)</b> |   |
|--|---|
| HCLP-27100-TS                              | High Compression Locking Peg, 2.7mm x 10mm        |
| HCLP-27120-TS                              | High Compression Locking Peg, 2.7mm x 12mm        |
| HCLP-27140-TS                              | High Compression Locking Peg, 2.7mm x 14mm        |
| HCLP-27160-TS                              | High Compression Locking Peg, 2.7mm x 16mm        |
| HCLP-27180-TS                              | High Compression Locking Peg, 2.7mm x 18mm        |
| HCLP-27190-TS                              | High Compression Locking Peg, 2.7mm x 19mm        |
| HCLP-27200-TS                              | High Compression Locking Peg, 2.7mm x 20mm        |
| HCLP-27210-TS                              | High Compression Locking Peg, 2.7mm x 21mm        |
| HCLP-27220-TS                              | High Compression Locking Peg, 2.7mm x 22mm        |
| HCLP-27230-TS                              | High Compression Locking Peg, 2.7mm x 23mm        |
| HCLP-27240-TS                              | High Compression Locking Peg, 2.7mm x 24mm        |
| HCLP-27260-TS                              | High Compression Locking Peg, 2.7mm x 26mm        |
| HCLP-27280-TS                              | High Compression Locking Peg, 2.7mm x 28mm        |
| HCLP-27300-TS                              | High Compression Locking Peg, 2.7mm x 30mm        |
| HCLP-27320-TS                              | High Compression Locking Peg, 2.7mm x 32mm        |
| <b>Threaded Pegs, Non-Locking (Ti)</b>     |   |
| TPNL-27100-TS                              | Threaded Peg, Non-Locking, 2.7mm x 10mm           |
| TPNL-27120-TS                              | Threaded Peg, Non-Locking, 2.7mm x 12mm           |
| TPNL-27140-TS                              | Threaded Peg, Non-Locking, 2.7mm x 14mm           |
| TPNL-27160-TS                              | Threaded Peg, Non-Locking, 2.7mm x 16mm           |
| TPNL-27180-TS                              | Threaded Peg, Non-Locking, 2.7mm x 18mm           |
| TPNL-27200-TS                              | Threaded Peg, Non-Locking, 2.7mm x 20mm           |
| TPNL-27220-TS                              | Threaded Peg, Non-Locking, 2.7mm x 22mm           |
| TPNL-27240-TS                              | Threaded Peg, Non-Locking, 2.7mm x 24mm           |
| TPNL-27260-TS                              | Threaded Peg, Non-Locking, 2.7mm x 26mm           |
| TPNL-27280-TS                              | Threaded Peg, Non-Locking, 2.7mm x 28mm           |
| TPNL-27300-TS                              | Threaded Peg, Non-Locking, 2.7mm x 30mm           |
| TPNL-27320-TS                              | Threaded Peg, Non-Locking, 2.7mm x 32mm           |
| <b>Polyaxial Screws, Locking (CoCr)</b>    |   |
| PALS-25100-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 10mm |
| PALS-25120-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 12mm |
| PALS-25140-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 14mm |
| PALS-25160-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 16mm |
| PALS-25180-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 18mm |
| PALS-25200-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 20mm |
| PALS-25220-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 22mm |
| PALS-25240-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 24mm |
| PALS-25260-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 26mm |
| PALS-25280-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 28mm |
| PALS-25300-CC                              | Screw, Cannulated Polyaxial Locking, 2.5mm x 30mm |
| <b>Cortical Screws, Non-Locking (Ti)</b>   |   |
| PANL-35080-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 8mm         |
| PANL-35090-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 9mm         |
| PANL-35100-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 10mm        |
| PANL-35110-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 11mm        |
| PANL-35120-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 12mm        |
| PANL-35130-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 13mm        |
| PANL-35140-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 14mm        |
| PANL-35150-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 15mm        |
| PANL-35160-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 16mm        |
| PANL-35180-TS                              | Screw, Cortical, Non-Locking, 3.5mm x 18mm        |

| <b>Cortical Screws, Locking (Ti)</b>           |   |
|--|---|
| COLS-35080-TS                                  | Screw, Cortical Locking, 3.5mm x 8mm                                  |
| COLS-35090-TS                                  | Screw, Cortical Locking, 3.5mm x 9mm                                  |
| COLS-35100-TS                                  | Screw, Cortical Locking, 3.5mm x 10mm                                 |
| COLS-35110-TS                                  | Screw, Cortical Locking, 3.5mm x 11mm                                 |
| COLS-35120-TS                                  | Screw, Cortical Locking, 3.5mm x 12mm                                 |
| COLS-35130-TS                                  | Screw, Cortical Locking, 3.5mm x 13mm                                 |
| COLS-35140-TS                                  | Screw, Cortical Locking, 3.5mm x 14mm                                 |
| COLS-35150-TS                                  | Screw, Cortical Locking, 3.5mm x 15mm                                 |
| COLS-35160-TS                                  | Screw, Cortical Locking, 3.5mm x 16mm                                 |
| COLS-35180-TS                                  | Screw, Cortical Locking, 3.5mm x 18mm                                 |
| <b>System Instrumentation</b>                  |   |
| DPGA-SMS-030<br>or<br>DPGA-UNV-030             | Depth Gauge, Sm. Standard, 30mm<br>or<br>Depth Gauge, Universal, 30mm |
| DRLL-SSC-20040                                 | Drill, Solid Side Cutting, 2.0mm x 40mm                               |
| DRLL-SSC-25040                                 | Drill, Solid Side Cutting, 2.5mm x 40mm                               |
| DRLL-PLS-20                                    | Drill, Cannulated Polyaxial Locking Screw, 2.0mm x 40mm               |
| DRVR-AOS-S20                                   | Driver, Peg, Torque Limiting  |
| DRVR-UQC-T10                                   | Driver, Universal Quick Connect, T10                                  |
| DRVR-AOS-PLS                                   | Driver, AO Connection, Polyaxial Locking Screw                        |
| HNDL-UQC-FXD<br>or<br>HNDL-AQC-FXD             | Handle, Universal Quick Connect, Fixed<br>or<br>Handle, AO QC, Fixed  |
| HNDL-SQC-FXD                                   | Handle, Small Quick Connect, Fixed                                    |
| FRCP-BHM-RTC                                   | Forceps, Bone Holding Medium, Ratcheting                              |
| TPDG-THD-DG20                                  | Thread-in Drill Guide, 2.0mm  |
| TPDG-THD-DG25                                  | Thread-in Drill Guide, 2.5mm  |
| TPDG-DSD-2025                                  | Tissue Protector / Drill Guide, Dual Sided, 2.0mm x 2.5mm             |
| GMN-ID-PLS                                     | Initial Driver, Polyaxial Locking Screw                               |
| GMN-CDG-PLS                                    | Cannulated Depth Gauge, Polyaxial Locking Screw                       |
| GMN-FSP-PLB                                    | GEMINUS Plate Bender  |
| GMN-FSP-PLH                                    | GEMINUS Plate Holder  |
| PDG-AIM-015                                    | AIMing Guides, 1.5mm  |
| PLS-AIM-0910                                   | PLS AIMing Guide  |
| KWIR-STD-09152                                 | K-wire, Standard Tip, .9mm x 152mm                                    |
| KWIR-STD-15127                                 | K-wire, Standard Tip, 1.5mm x 127mm                                   |
| GMN-HP-DG15                                    | GEMINUS Hook Plate, Reduction Tool                                    |
| <b>Sterilization Trays</b>                     |   |
| GMN-FSP-TRAY                                   | GEMINUS Volar Distal Radius Plate System Sterilization Tray           |
| GMN-ACC-MOD1                                   | GEMINUS Volar Distal Radius Plate System Accessory Module 1           |
| GMN-TRAY-SST                                   | GEMINUS Volar Distal Radius Plate System Sterilization Tray, 304      |
| GMN-MOD1-SST                                   | GEMINUS Volar Distal Radius Plate System Accessory Module 1, 304      |
| <b>Optional Drill Block System (U.S. ONLY)</b> |   |
| GMN-DBK-RTS                                    | GEMINUS Drill Block, Right, Standard                                  |
| GMN-DBK-LTS                                    | GEMINUS Drill Block, Left, Standard                                   |
| GMN-DBK-RTW                                    | GEMINUS Drill Block, Right, Wide                                      |
| GMN-DBK-LTW                                    | GEMINUS Drill Block, Left, Wide                                       |
| GMN-DBK-RTN                                    | GEMINUS Drill Block, Right, Narrow                                    |
| GMN-DBK-LTN                                    | GEMINUS Drill Block, Left, Narrow                                     |
| TPDG-DBK-DG20                                  | GEMINUS Drill Block, Drill Guide, 2.0mm                               |
| DBK-AIM-015                                    | GEMINUS Drill Block, AIMing Guide, 1.5mm                              |
| GMN-DBK-MOD                                    | GEMINUS Drill Block, Sterilization Module                             |

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