

## **GEMINUS<sup>®</sup> Distal Radius Plating System**

### **INSTRUCTIONS FOR USE FOR EAST & SOUTHEAST ASIA**

---

**R:** For use by physicians only. Federal Law (USA) restricts this device to sale by or on order of a physician.

#### **Failure to follow instructions may lead to patient injury.**

This package insert is designed to provide Instructions for Use of the GEMINUS<sup>®</sup> Distal Radius Plating System.

#### **Description**

The Skeletal Dynamics GEMINUS<sup>®</sup> Distal Radius Plating System contains bone plates for the repair of distal volar radial fractures. Included in the set are titanium bone screws, fixation pegs, fragment plates, and specialized instrumentation. Also included are a Hook Plate Extension to buttress a volar marginal fragment, and cannulated cobalt chrome polyaxial locking screws for trajectories different than those of the fixed angled bone plates.

The GEMINUS<sup>®</sup> Distal Radius Plates are available in various sizes and are made of medical grade titanium alloy. Cortical screws affix the plate to the diaphysis and fixed angle pegs are used for distal bone fragments. The system is provided non-sterile and is sterilized in the user facility.

The GEMINUS<sup>®</sup> Distal Radius Plating System is comprised of:

- Titanium alloy plates, washers and screws
- CoCr Cannulated Polyaxial Locking Screw (PLS)
- Stainless steel K-wires (for provisional fixation; not for implantation)
- System specific instrumentation

#### **Indications**

The GEMINUS<sup>®</sup> Distal Radius Plating System is intended for the fixation of fractures and osteotomies of the distal radius.

#### **Contraindications**

Prior to using the GEMINUS<sup>®</sup> Distal Radius Plating System, ensure that none of the following patient conditions are present: active or latent infection, insufficient quantity or quality of bone and/or soft tissue, material sensitivity, or patients who are unwilling or incapable of following postoperative care instructions.

#### **⚠ Warnings and Precautions**

- All screws must be implanted and fully tightened into the plate to maintain the integrity and strength of the finished construct. If the screws are not attached and/or fully tightened, a non-union, delayed union or construct failure may occur.
- The use of power tools for the installation of the screws and pegs is not recommended and may lead to cross threading and damage to the screws and/or plates.
- The information in this document should be shared with the patient.
- The patient should be informed about the importance of following the post operative rehabilitation prescribed to fully understand the limitations in activities of daily living. The patient must be warned that failure to follow postoperative care instructions may cause the implant or treatment to fail.

- Potential GEMINUS<sup>®</sup> Distal Radius Plating System construct failures such as implant breakage, loosening of the construct and/or fixation, delayed union, or non-union may occur as a result of non-compliance to postoperative rehabilitation, excessive wrist activities or construct overloading.
- **DO NOT** reuse any of the GEMINUS<sup>®</sup> Distal Radius Plating System implantable components. Reuse may compromise the structural integrity of the construct and/or lead to failure or infection, which may result in patient injury.
- **DO NOT** open the volar capsule as it may devascularize fracture fragments and destabilize the volar wrist ligaments.
- Use only one 2.7mm Peg (High Compression or Fully Threaded, Non-Locking) in each head of the GEMINUS<sup>®</sup> Distal Radius Plate.
- Use only one 2.5mm PLS in each head of the GEMINUS<sup>®</sup> Distal Radius Plate.
- **DO NOT** use the PLS in the most distal hole(s) on the lunate head of the GEMINUS<sup>®</sup> Distal Radius Plate.
- GEMINUS<sup>®</sup> Drill Blocks are only compatible with GEMINUS<sup>®</sup> Distal Radius Plates containing a Gold PDG in the shaft.
- Protect the GEMINUS<sup>®</sup> Distal Radius Plating System's implantable components against scratching or nicking. Such stress concentration can lead to implant failure.
- Before using the GEMINUS<sup>®</sup> Distal Radius Plating System, inspect all implants and instruments for wear, disfiguration and physical damage. If evidence of wear, disfiguration or physical damage is found, **DO NOT** use and contact your local Skeletal Dynamics representative or the Skeletal Dynamics Customer Care Department.
- Assure Peg Driver tip does not show any signs of wear or distress such as rounded square edges, excessive depth marks from peg recess insertion, or deformed twisted tip. If such evidence is found for Peg Driver, **DO NOT USE** and contact your local Skeletal Dynamics representative or the Skeletal Dynamics Customer Care Department for replacement.
- **DO NOT** permanently implant the Skeletal Dynamics K-Wires; they are intended to be used during provisional fixation of the GEMINUS<sup>®</sup> Distal Radius Plate.
- **DO NOT** permanently implant the pre-loaded Drill Guides, Drill Blocks, or AIMing Guides; they are intended to be removed prior to peg insertion.
- **DO NOT** use peg/screw lengths that will excessively protrude through the far cortex as it may result in soft tissue irritation.
- The maximum angulation of the PLS should not exceed 10° from the trajectory of the respective hole.
- The Non-Locking Threaded Pegs are NOT intended to provide subchondral support. Their use should be limited to capture remote bone fragments where partially or fully threaded pegs cannot be used.
- The Skeletal Dynamics GEMINUS<sup>®</sup> Distal Radius Plating System is to be used only with Skeletal Dynamics instruments, implants and accessories.
- Dispose of contaminated implants and instruments per established facility guidelines and protocols.
- Caution should be taken for interference to pacemakers during electrocautery or by uncertified drill power sources.
- Seek medical help immediately if implant malfunctions.

- To maintain traceability of the GEMINUS<sup>®</sup> Distal Radius Plating System implantable components, you must record each of the respective components LOT numbers into the patient medical records post implantation.

### **Potential Adverse Events**

The following are potential risks that have been associated with wrist surgery: infection, nonunion, nerve or soft tissue damage, persistent pain, stiffness of the fingers, loosening or migration of the implants resulting in misalignment.

### **⚠ MR Safety Information.**

The GEMINUS<sup>®</sup> Distal Radius Plating System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the GEMINUS<sup>®</sup> Distal Radius Plating System in the MR environment is unknown. Scanning a patient who has this device implanted may result in patient injury.

### **Directions for Use**

The GEMINUS<sup>®</sup> Distal Radius Plating System should only be used by surgeons who have experience with this system. Each surgeon must evaluate the appropriateness for the use of the GEMINUS<sup>®</sup> Distal Radius Plating System based on their clinical experiences.

Please refer to the *GEMINUS<sup>®</sup> Distal Radius Plating System's Surgical Technique Guide* (MKT-00005-00) to review the surgical approach as described by Jorge L. Orbay, M.D. of the *Miami Hand and Upper Extremity Institute* located in Miami, Florida, USA.

### **Cleaning, Sterilization, and Inspection**

For instructions on cleaning, disinfection, sterilization and inspection of the Geminus<sup>®</sup> products please refer to cleaning and sterilization instructions for use (IFU-04056-00). Refer to the *Geminus Distal Radius System Surgical Technique Guide* (MKT-00005-00) for proper kitted tray arrangement.

### **Resources**

The latest version of the Instructions for Use may be requested in a physical format by email ([orders@skeletaldynamics.com](mailto:orders@skeletaldynamics.com)) or by phone (+1-877-753-5396). The physical copy will be provided within 7 calendar days of receiving a request from the user or at the time of delivery of the device if so, requested at the time of order.

For the most current instructions for use visit [www.skeletaldynamics.com/resources](http://www.skeletaldynamics.com/resources). Instructions for Use should always be reviewed before using or implanting a device.

### **Disclaimer of Warranty and Limited Remedies**

Skeletal Dynamics, Inc. makes no express or implied warranty, including any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described in this publication. Skeletal Dynamics, Inc. shall not be liable under any circumstances for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has authority to bind Skeletal Dynamics, Inc. to any representation or warranty except as specifically set forth in this publication. Descriptions or specifications provided by Skeletal Dynamics, Inc. in any publication are only included to generally describe the product when manufactured and do not constitute any express warranties



Skeletal Dynamics, Inc  
7300 N. Kendall Dr. / Suite 800  
Miami, FL 33156  
United States of America  
1-877-753-5396

**LRP:**

**Hong Kong Representative: Emergo  
Hong Kong Limited**  
18/F, Delta House  
3 ON YIU Street  
Shatin, NT HONG KONG

Contact No: +852 2276 9196

HKMD No. 251722  
没有中文版使用说明



**Malaysia Representative:**

**Emergo Malaysia Sdn. Bhd.**









Level 16, 1 Sentral  
Jalan Stesen Sentral  
5 KL Sentral, 50470  
Kuala Lumpur Malaysia  
Website: [www.emergobyul.com](http://www.emergobyul.com)




MDA Registration #:

GC7785025-201325

# SYMBOLS GLOSSARY

Symbols that follow BS EN ISO 15223-1 / Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied.

Symbol	Symbol Reference Number and Title	Description
	5.2.7 Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.
	5.1.6 Catalogue Number	Indicates the manufacturer's catalogue number so the medical device can be identified.
	5.4.3 Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.
	5.1.5 Batch Code	Indicates the manufacturer's batch code or lot can be identified.
	5.4.2 Do Not Re-use	Indicates a medical device that is intended for one single use only.
	5.1.3 Date of Manufacture	Indicates the date when the medical device was manufactured.
	5.1.1 Manufacturer	Indicates the medical device manufacturer.
	5.2.8 Do Not Use if Package is Damaged	Indicates a medical device that should not be used if the package has been damaged or opened.

Symbol	Symbol Reference Number and Title	Description
	<p>5.1.4 Use-by-date</p>	<p>Indicates the date after which the medical device is not to be used.</p>
	<p>5.7.10 Unique Device Identifier</p>	<p>Indicates a carrier that contains unique device identifier information.</p>
	<p>5.7.7 Medical device</p>	<p>Indicated the item is a medical device</p>

## GEMINUS® Distal Radius Plating System Inventory Control Sheet

GEMINUS Distal Radius Distal Radius Plates (Ti)			
GEMINUS Plate, Narrow, 3 Hole, Left GMN-LTN-3HL (01)00841506101569 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101569</small> </div>	GEMINUS Plate, Narrow, 4 Hole, Left GMN-LTN-4HL (01)00841506101576 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101576</small> </div>		
GEMINUS Plate, Narrow, 3 Hole, Right GMN-RTN-3HL (01)00841506101620 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101620</small> </div>	GEMINUS Plate, Narrow, 4 Hole, Right GMN-RTN-4HL (01)00841506101637 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101637</small> </div>		
GEMINUS Plate, Standard, 3 Hole, Left GMN-LTS-3HL (01)00841506101583 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101583</small> </div>	GEMINUS Plate, Standard, 4 Hole, Left GMN-LTS-4HL (01)00841506101590 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101590</small> </div>		
GEMINUS Plate, Standard, 3 Hole, Right GMN-RTS-3HL (01)00841506101644 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101644</small> </div>	GEMINUS Plate, Standard, 4 Hole, Right GMN-RTS-4HL (01)00841506101651 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101651</small> </div>		
GEMINUS Plate, Standard, 7 Hole, Left GMN-LTS-7HL (01)00841506101606 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101606</small> </div>	GEMINUS Plate, Wide, 4 Hole, Left GMN-LTW-4HL (01)00841506101613 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101613</small> </div>		
GEMINUS Plate, Standard, 7 Hole, Right GMN-RTS-7HL (01)00841506101668 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101668</small> </div>	GEMINUS Plate, Wide, 4 Hole, Right GMN-RTW-4HL (01)00841506101675 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101675</small> </div>		
PROTEAN Fragment Plates (Ti)			
PROTEAN Fragment Plate, Double Hockey Stick PRT-FSP-LR (01)00841506102917 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506102917</small> </div>	PROTEAN Fragment Plate, Y PRT-FSP-YS (01)0084150610291 <div style="text-align: right; margin-top: 10px;">   <small>(01) 0084150610291</small> </div>		
PROTEAN Fragment Plate, Distal Ulna PRT-FSP-DU (01)00841506102900 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506102900</small> </div>	PROTEAN Fragment Plate, Radial Column Plate, Left PRT-RCP-LT (01)00841506109930 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506109930</small> </div>		
PROTEAN Fragment Plate, Radial Column Plate, Right PRT-RCP-RT (01)00841506109923 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506109923</small> </div>	PROTEAN Fragment Plate, Central Column Plate, Left PRT-CCP-LT (01)00841506109954 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506109954</small> </div>		
PROTEAN Fragment Plate, Central Column Plate, Right PRT-CCP-RT (01)00841506109947 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506109947</small> </div>			
GEMINUS Hook Plate (Ti)			
GEMINUS Hook Plate GMN-HP (01)00841506101514 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101514</small> </div>	GEMINUS Hook Plate Screw GMN-HP-SCRW (01)0084150610154 <div style="text-align: right; margin-top: 10px;">   <small>(01) 00841506101545</small> </div>		









### Single Use (Disposable) Instruments







AIMing Guides, 1.5mm PDG-AIM-015 (01)00841506102870	PLS AIMing Guides, .9mm x 10° (Cone Tool) PLS-AIM-0910 (01)00841506102887
 <small>(01) 00841506102870</small>	 <small>(01) 00841506102887</small>
K-Wire, Standard Tip, 0.9mm x 152mm KWIR-STD-09152 (01)00841506102498	K-Wire, Standard Tip, 1.6mm x 127mm KWIR-STD-15127 (01)00841506102504
 <small>(01) 00841506102498</small>	 <small>(01) 00841506102504</small>

### Smooth Pegs, Locking (Ti) - Gray

Smooth Peg, Locking, 2.0mm x 10mm SPLS-20100-TS (01)00841506102948	Smooth Peg, Locking, 2.0mm x 20mm SPLS-20200-TS (01)00841506103013
 <small>(01) 00841506102948</small>	 <small>(01) 00841506103013</small>
Smooth Peg, Locking, 2.0mm x 12mm SPLS-20120-TS (01)00841506102955	Smooth Peg, Locking, 2.0mm x 21mm SPLS-20210-TS (01)00841506103020
 <small>(01) 00841506102955</small>	 <small>(01) 00841506103020</small>
Smooth Peg, Locking, 2.0mm x 14mm SPLS-20140-TS (01)00841506102962	Smooth Peg, Locking, 2.0mm x 22mm SPLS-20220-TS (01)00841506103037
 <small>(01) 00841506102962</small>	 <small>(01) 00841506103037</small>
Smooth Peg, Locking, 2.0mm x 16mm SPLS-20160-TS (01)00841506102979	Smooth Peg, Locking, 2.0mm x 23mm SPLS-20230-TS (01)00841506103044
 <small>(01) 00841506102979</small>	 <small>(01) 00841506103044</small>
Smooth Peg, Locking, 2.0mm x 17mm SPLS-20170-TS (01)00841506102986	Smooth Peg, Locking, 2.0mm x 24mm SPLS-20240-TS (01)00841506103051
 <small>(01) 00841506102986</small>	 <small>(01) 00841506103051</small>
Smooth Peg, Locking, 2.0mm x 18mm SPLS-20180-TS (01)00841506102993	Smooth Peg, Locking, 2.0mm x 26mm SPLS-20260-TS (01)00841506103068
 <small>(01) 00841506102993</small>	 <small>(01) 00841506103068</small>
Smooth Peg, Locking, 2.0mm x 19mm SPLS-20190-TS (01)00841506103006	Smooth Peg, Locking, 2.0mm x 28mm SPLS-20280-TS (01)00841506103075
 <small>(01) 00841506103006</small>	 <small>(01) 00841506103075</small>

### High Compression Locking Pegs (Ti) - Bronze





High Compression Locking Peg, 2.7mm x 10mm HCLP-27100-TS (01)00841506101682	High Compression Locking Peg, 2.7mm x 21mm HCLP-27210-TS (01)00841506101750
 <small>(01) 00841506101682</small>	 <small>(01) 00841506101750</small>
High Compression Locking Peg, 2.7mm x 12mm HCLP-27120-TS (01)00841506101699	High Compression Locking Peg, 2.7mm x 22mm HCLP-27220-TS (01)00841506101767
 <small>(01) 00841506101699</small>	 <small>(01) 00841506101767</small>
High Compression Locking Peg, 2.7mm x 14mm HCLP-27140-TS (01)00841506101705	High Compression Locking Peg, 2.7mm x 23mm HCLP-27230-TS (01)00841506101774
 <small>(01) 00841506101705</small>	 <small>(01) 00841506101774</small>
High Compression Locking Peg, 2.7mm x 16mm HCLP-27160-TS (01)00841506101712	High Compression Locking Peg, 2.7mm x 24mm HCLP-27240-TS (01)00841506101781
 <small>(01) 00841506101712</small>	 <small>(01) 00841506101781</small>

High Compression Locking Peg, 2.7mm x 18mm HCLP-27180-TS (01)00841506101729	 (01) 00841506101729	High Compression Locking Peg, 2.7mm x 26mm HCLP-27260-TS (01)00841506101798	 (01) 00841506101798
High Compression Locking Peg, 2.7mm x 19mm HCLP-27190-TS (01)00841506101736	 (01) 00841506101736	High Compression Locking Peg, 2.7mm x 28mm HCLP-27280-TS (01)00841506101804	 (01) 00841506101804
High Compression Locking Peg, 2.7mm x 20mm HCLP-27200-TS (01)00841506101743	 (01) 00841506101743	High Compression Locking Peg, 2.7mm x 30mm HCLP-27300-TS (01)00841506101811	 (01) 00841506101811















**Threaded Pegs, Locking (Ti) - Pink**

Threaded Peg, Locking, 2.3mm x 10mm TPLS-23100-TS (01)00841506103358	 (01) 00841506103358	Threaded Peg, Locking, 2.3mm x 21mm TPLS-23210-TS (01)00841506103433	 (01) 00841506103433
Threaded Peg, Locking, 2.3mm x 12mm TPLS-23120-TS (01)00841506103365	 (01) 00841506103365	Threaded Peg, Locking, 2.3mm x 22mm TPLS-23220-TS (01)00841506103440	 (01) 00841506103440
Threaded Peg, Locking, 2.3mm x 14mm TPLS-23140-TS (01)00841506103372	 (01) 00841506103372	Threaded Peg, Locking, 2.3mm x 23mm TPLS-23230-TS (01)00841506103457	 (01) 00841506103457
Threaded Peg, Locking, 2.3mm x 16mm TPLS-23160-TS (01)00841506103389	 (01) 00841506103389	Threaded Peg, Locking, 2.3mm x 24mm TPLS-23240-TS (01)00841506103464	 (01) 00841506103464
Threaded Peg, Locking, 2.3mm x 17mm TPLS-23170-TS (01)00841506103396	 (01) 00841506103396	Threaded Peg, Locking, 2.3mm x 26mm TPLS-23260-TS (01)00841506103471	 (01) 00841506103471
Threaded Peg, Locking, 2.3mm x 18mm TPLS-23180-TS (01)00841506103402	 (01) 00841506103402	Threaded Peg, Locking, 2.3mm x 28mm TPLS-23280-TS (01)00841506103488	 (01) 00841506103488
Threaded Peg, Locking, 2.3mm x 19mm TPLS-23190-TS (01)00841506103419	 (01) 00841506103419	Threaded Peg, Locking, 2.3mm x 30mm TPLS-23300-TS (01)00841506103495	 (01) 00841506103495
Threaded Peg, Locking, 2.3mm x 20mm TPLS-23200-TS (01)00841506103426	 (01) 00841506103426		

**Threaded Pegs, Non-Locking (Ti) - Blue**

Threaded Peg, Non-Locking, 2.7mm x 10mm TPNL-27100-TS (01)00841506103518	 (01) 00841506103518	Threaded Peg, Non-Locking, 2.7mm x 22mm TPNL-27220-TS (01)00841506103570	 (01) 00841506103570
Threaded Peg, Non-Locking, 2.7mm x 12mm TPNL-27120-TS (01)00841506103525	 (01) 00841506103525	Threaded Peg, Non-Locking, 2.7mm x 24mm TPNL-27240-TS (01)00841506103587	 (01) 00841506103587

Threaded Peg, Non-Locking, 2.7mm x 14mm TPNL-27140-TS (01)00841506103532  (01) 00841506103532	Threaded Peg, Non-Locking, 2.7mm x 26mm TPNL-27260-TS (01)00841506103594  (01) 00841506103594
Threaded Peg, Non-Locking, 2.7mm x 16mm TPNL-27160-TS (01)00841506103549  (01) 00841506103549	Threaded Peg, Non-Locking, 2.7mm x 28mm TPNL-27280-TS (01)00841506103600  (01) 00841506103600
Threaded Peg, Non-Locking, 2.7mm x 18mm TPNL-27180-TS (01)00841506103556  (01) 00841506103556	Threaded Peg, Non-Locking, 2.7mm x 30mm TPNL-27300-TS (01)00841506103617  (01) 00841506103617
Threaded Peg, Non-Locking, 2.7mm x 20mm TPNL-27200-TS (01)00841506103563  (01) 00841506103563	
<b>Cannulated Polyaxial Screws, Locking (CoCr) - Chrome</b>	
Screw, Polyaxial Locking, 2.5mm x 10mm, Cannulated PALS-25100-CC (01)00841506102665  (01) 00841506102665	Screw, Polyaxial Locking, 2.5mm x 22mm, Cannulated PALS-25220-CC (01)00841506102726  (01) 00841506102726
Screw, Polyaxial Locking, 2.5mm x 12mm, Cannulated PALS-25120-CC (01)00841506102672  (01) 00841506102672	Screw, Polyaxial Locking, 2.5mm x 24mm, Cannulated PALS-25240-CC (01)00841506102733  (01) 00841506102733
Screw, Polyaxial Locking, 2.5mm x 14mm, Cannulated PALS-25140-CC (01)00841506102689  (01) 00841506102689	Screw, Polyaxial Locking, 2.5mm x 26mm, Cannulated PALS-25260-CC (01)00841506102740  (01) 00841506102740
Screw, Polyaxial Locking, 2.5mm x 16mm, Cannulated PALS-25160-CC (01)00841506102696  (01) 00841506102696	Screw, Polyaxial Locking, 2.5mm x 28mm, Cannulated PALS-25280-CC (01)00841506102757  (01) 00841506102757
Screw, Polyaxial Locking, 2.5mm x 18mm, Cannulated PALS-25180-CC (01)00841506102702  (01) 00841506102702	Screw, Polyaxial Locking, 2.5mm x 30mm, Cannulated PALS-25300-CC (01)00841506102764  (01) 00841506102764
Screw, Polyaxial Locking, 2.5mm x 20mm, Cannulated PALS-25200-CC (01)00841506102719  (01) 00841506102719	
<b>Cortical Screws, Locking (Ti) - Bronze</b>	
Screw, Cortical Locking, 3.5mm x 8mm COLS-35080-TS (01)00841506101071  (01) 00841506101071	Screw, Cortical Locking, 3.5mm x 13mm COLS-35130-TS (01)00841506101125  (01) 00841506101125
Screw, Cortical Locking, 3.5mm x 9mm COLS-35090-TS (01)00841506101088  (01) 00841506101088	Screw, Cortical Locking, 3.5mm x 14mm COLS-35140-TS (01)00841506101132  (01) 00841506101132
Screw, Cortical Locking, 3.5mm x 10mm COLS-35100-TS (01)00841506101095  (01) 00841506101095	Screw, Cortical Locking, 3.5mm x 15mm COLS-35150-TS (01)00841506101149  (01) 00841506101149

Screw, Cortical Locking, 3.5mm x 11mm COLS-35110-TS (01)00841506101101  (01)00841506101101	Screw, Cortical Locking, 3.5mm x 16mm COLS-35160-TS (01)00841506101156  (01)00841506101156
Screw, Cortical Locking, 3.5mm x 12mm COLS-35120-TS (01)00841506101118  (01)00841506101118	Screw, Cortical Locking, 3.5mm x 18mm COLS-35180-TS (01)00841506101163  (01)00841506101163
<b>Cortical Screws, Non-Locking (Ti) - Blue</b>	
Screw, Cortical Non-Locking, 3.5mm x 8mm PANL-35080-TS (01)00841506102771  (01)00841506102771	Screw, Cortical Non-Locking, 3.5mm x 13mm PANL-35130-TS (01)00841506102825  (01)00841506102825
Screw, Cortical Non-Locking, 3.5mm x 9mm PANL-35090-TS (01)00841506102788  (01)00841506102788	Screw, Cortical Non-Locking, 3.5mm x 14mm PANL-35140-TS (01)00841506102832  (01)00841506102832
Screw, Cortical Non-Locking, 3.5mm x 10mm PANL-35100-TS (01)00841506102795  (01)00841506102795	Screw, Cortical Non-Locking, 3.5mm x 15mm PANL-35150-TS (01)00841506102849  (01)00841506102849
Screw, Cortical Non-Locking, 3.5mm x 11mm PANL-35110-TS (01)00841506102801  (01)00841506102801	Screw, Cortical Non-Locking, 3.5mm x 16mm PANL-35160-TS (01)00841506102856  (01)00841506102856
Screw, Cortical Non-Locking, 3.5mm x 12mm PANL-35120-TS (01)00841506102818  (01)00841506102818	Screw, Cortical Non-Locking, 3.5mm x 18mm PANL-35180-TS (01)00841506102863  (01)00841506102863
<b>Washers (Ti) - Blue</b>	
Washer, Inside Ø2.7mm, Outside Ø5.0mm WBTN-2750-T (01)00841506103730  (01)00841506103730	

## Geminus® Distal Radius System Instruments

Catalogue	Description
DBK-AIM-015	GEMINUS Drill Block, AIMing Guide, 1.5mm
DPGA-SMS-030	Depth Gauge, Sm. Standard, 30mm
DPGA-UNV-030	Depth Gauge, Universal, 30mm
DRLL-20040R	Drill, Solid Side Cutting, 2.0mm x 40mm, Reusable
DRLL-25040R	Drill, Solid Side Cutting, 2.5mm x 40mm, Reusable
DRLL-PLS-20	Drill, Cannulated, Polyaxial Locking Screw, 2.0mm
DRLL-PLS-20R	Drill, Cannulated, Polyaxial Locking Screw, 2.0mm Reusable
DRLL-SSC-20040	Drill, Solid Side Cutting, 2.0mm x 40mm
DRLL-SSC-25040	Drill, Solid Side Cutting, 2.5mm x 40mm
DRVR-AOS-PLS	Driver, AO Connection, Polyaxial Locking Screw
DRVR-AOS-PLSR	Driver, AO Connection, Polyaxial Locking Screw, Reusable
DRVR-AOS-S20	Driver, Peg, Torque Limiting
DRVR-AOS-S20R	Driver, Peg, Torque Limiting, Reusable
DRVR-S20	Driver, Peg
DRVR-S20R	Driver, Peg, Reusable
DRVR-T10R	Driver, Universal QC, T10, Reusable
DRVR-UQC-T10	Driver, Universal Quick Connect, T10
FRCP-BHM-RTC	Forceps, Bone Holding Medium, Ratcheting
GMN-ACC-MOD1	GEMINUS Volar Distal Radius Plate System Accessory Module 1
GMN-CDG-PLS	Cannulated Depth Gauge, Polyaxial Locking Screw
GMN-DBK-LTN	GEMINUS Drill Block, Left, Narrow
GMN-DBK-LTS	GEMINUS Drill Block, Left, Standard
GMN-DBK-LTW	GEMINUS Drill Block, Left, Wide
GMN-DBK-MOD	GEMINUS Volar Distal Radius Plate System Drill Block Module
GMN-DBK-RTN	GEMINUS Drill Block, Right, Narrow
GMN-DBK-RTS	GEMINUS Drill Block, Right, Standard
GMN-DBK-RTW	GEMINUS Drill Block, Right, Wide
GMN-DBK-SYS	GEMINUS Drill Block System
GMN-FSP-PLB	GEMINUS Plate Bender
GMN-FSP-PLH	GEMINUS Plate Holder
GMN-FSP-TRAY	GEMINUS Volar Distal Radius Plate System Sterilization Tray
GMN-HP-DG15	GEMINUS Hook Plate, Reduction Tool
GMN-ID-PLS	Initial Driver, Polyaxial Locking Screw
GMN-LNG-MOD	GEMINUS Volar Distal Radius Plate System Long Plate Module
GMN-TRAY-BASE	GEMINUS Volar Distal Radius Plate System Base
GMN-TRAY-CAD	GEMINUS Volar Distal Radius Plate System Screw Caddy
GMN-TRAY-INSRT	GEMINUS Volar Distal Radius Plate System Instrument Tray Insert
GMN-TRAY-LID	GEMINUS Volar Distal Radius Plate System Lid
HNDL-AQC-FXD	Handle, AO QC, Fixed
HNDL-PDL-FX1	Handle, Paddle, Fixed
HNDL-PDL-RT1	Handle, Paddle, Ratcheting
HNDL-SQC-FX1	Handle, Small QC, Fixed
HNDL-SQC-FXD	Handle, Small QC, Fixed
HNDL-SQC-RAT	Handle, Small QC, Ratcheting
HNDL-UQC-FXD	Handle, Universal Quick Connect, Fixed
INST-KPE-STD	Instrument, Key Periosteal Elevator, Standard

INST-MHR-STD	Instrument, Mini-Hohmann Retractor, Standard
INST-SHP-STD	Instrument, Sharp Hook Probe, Standard
KWIR-STD-09152	K-Wire, Standard Tip, .9mm x 152mm
KWIR-STD-15127	K-Wire, Standard Tip, 1.6mm x 127mm
PDG-AIM-015	AIMing Guides, 1.5mm
PLS-AIM-0910	PLS AIMing Guides
PRT-BND-PLR	Protean Plate Bending Pliers
TPDG-DSD-2025	Tissue Protector / Drill Guide, Dual Sided, 2.0mm x 2.5mm
TPDG-GTS-25	Tissue Protector/Drill Guide, Slot, 2.5mm
TPDG-THD-DG20	Thread-in Drill Guide, 2.0mm
TPDG-THD-DG25	Thread-in Drill Guide, 2.5mm
TQH-S075	Torque Limiting Handle with Spin Cap, Mini, Jeweler