


IMPLATE[®]

wrist arthodesis nail

by skeletal dynamics

The First
Intramedullary Nail
System for Total
Wrist Fusion.



+ Less Invasive

+ Strong and Stable Construct

+ Precise Anatomic Alignment

IMPLATE[®]

wrist arthodesis nail




by skeletal dynamics

Unicortical Locking Screws

+ for nail stability

Threaded Slots

+ Provide up to 8mm of translation for compression or distraction




MicroSplines

+ Allow for rotational adjustments
+ Provide construct stability

Positional Adjustments

+ Rotation
+ Flexion-Extension
+ Radio-Ulnar Deviation



STRONG AND STABLE CONSTRUCT