


IJS[®]-ELBOW

elbow stabilization system
by skeletal dynamics

Comprehensive Elbow Stabilization



-
- + Stabilizes the unstable elbow
 - + Permits early active motion
 - + Reduces inconvenience and complications of external fixation devices
-

IJS[®]-ELBOW

elbow stabilization system
by skeletal dynamics

Proximal Locking Joint

Proximal Locking Screw

Distal Locking Joint

Distal Locking Screw

Sliding Slot

Compression Screw
(Polyaxial Non Locking)

Proximal
Connecting
Rod

Distal
Connecting
Rod

Axis Pin

Base Plate