



PROTEAN[®]

radial head plate
by skeletal dynamics


Minimal Footprint
Maximum Fixation



PROTEAN[®]

radial head plate

by skeletal dynamics



+ In-situ plate contouring

+ Various screw options for fixed angle or multi-directional applications

+ Malleable in all 3 planes




Plate Shaping


Vertical Plane:
Bending up to 30°
between any two nodes.



Horizontal Plane:
bending up to 5°
between any two nodes.



Transverse Plane:
bending up to 45°
between any two nodes.



Cutting to length:
Bend past 120° between
two adjacent nodes.

